Généralités

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© 2019 Hamilton Bonaduz AG. All rights reserved.

REF 238257/04 — 08/2019

HAMILTON®

Web: www.hamiltoncompany.com

USA: **800-648-5950**

Europe: **+41-58-610-10-10**


Hamilton Americas & Pacific Rim
4970 Energy Way
Reno, Nevada 89502 USA
Tel: +1-775-858-3000
Fax: +1-775-856-7259
sales@hamiltoncompany.com

Hamilton Europe, Asia & Africa
Via Crusch 8
CH-7402 Bonaduz, Switzerland
Tel: +41-58-610-10-10
Fax: +41-58-610-00-10
contact.pa.ch@hamilton.ch

To find a representative in your area, please visit hamiltoncompany.com

This guide may be available in other languages.
Visit www.hamiltoncompany.com for more information.

Operating Instructions for pH and ORP electrodes

Operating Instruction for pH and ORP electrodes

Introduction

These operating instructions are intended for all Hamilton pH and ORP sensors. The specifications must not be exceeded.

These instructions should be read, understood and followed by all staff using the device. Hamilton can assume no responsibility for damage and operational disruptions arising from failure to observe these instructions.

Liability

The Liability of Hamilton Bonaduz AG is detailed in the document «General Terms and Conditions of Sale and Delivery (GTS)», chapter 12.

Hamilton is expressly not liable for direct or indirect losses arising from use of the sensors. It must in particular be insured in this conjunction that malfunctions can occur on account of the inherently limited useful life of electrodes contingent upon their relevant applications. The user is responsible for the calibration, maintenance and regular replacement of the electrodes. In the case of critical electrode applications, Hamilton recommends using back-up measuring points in order to avoid consequential damages. The user is responsible for taking suitable precautions in the event of a electrode failure.

Intended Use

These electrode are intended for the measurement of pH or ORP values. For accurate measurements the surface of the sensor must be clean. A thorough rinsing with deionized water is recommended to remove any residue. pH and ORP sensors should only remain in the sample as long the measurement takes place. When the sensor is not in use it should be removed from the sample and stored in 3M KCl or Storage Solution.

Initial operation

Check the sensor for any damage at the time of un-packing. A damaged sensor should be returned to your Hamilton partner in the original packaging.

Every product for shipment or sent back for repair must be decontaminated. If working with hazardous liquids observe and carry out the maintenance procedures, paying particular attention to cleaning and decontamination. If the product becomes contaminated with biohazardous, radioactive or chemical material, it should be cleaned.

Preparing the Electrode

Carefully remove watering cap. Flush electrode with water.

Refillable electrodes: Prior to measuring, open stopper; close it after measuring. Check level of electrolyte in electrode. Minimum level is 20 mm below the refill opening. At the time of measurement the level of electrolyte must be higher than the level of the sample solution.

pH: Check inside of pH glass membrane for unwanted air bubbles. Shaking the electrode gently will cause any bubbles to rise to the top.

Sleeve diaphragm: Loosen the sleeve (plastic part) of the ground sleeve diaphragm by twisting upwards slightly. Allow a drop of electrolyte to flow out, then retighten the sleeve so that the electrolyte no longer flows out. Attention: Retighten sleeve with slight pressure only.

Should it not be possible to loosen the joint, dip the lower part of the electrode into water warmed to 40 - 50 °C for a few minutes, since the sleeve can be more easily loosened through warmth.

The loosening of the joint should be repeated occasionally to renew the electrolyte film in the diaphragm.

The flow of electrolyte increases significantly at higher

temperatures. This is caused by the varying expansion of the glass shaft and the plastic sleeve. Should the reservoir empty after a short period, or measurements have to be made at higher temperatures (50 °C) over a long period, tighten the sleeve after the working temperature has been reached.

Calibration and Measurement

For calibration, always use fresh Hamilton DuraCal pH or ORP buffers.

Air/gas bubbles may stick to the sensitive area of the sensor. As a consequence, the measurement value might be wrong. In case of a shutdown of the power supply (230V) the measurement value could be wrong.

Storing the Electrode

Electrodes should be stored with the watering cap on, containing 3M KCl solution (Ref 238036), Protelyte (Ref 238038) or Storage Solution (Ref 238931). Electrodes stored dry will exhibit unstable values for a short period. If the electrode inadvertently dries out, it can be regenerated by placing it overnight in Storage Solution.

Cleaning

In general, you can use acids, alkaline solutions and commonly used solvents for brief periods to clean electrodes with a glass shaft. Flush with water after cleaning. After cleaning, electrodes are likely to exhibit sluggish response times for a certain period, so place it in storage solution for 15 minutes after cleaning.

Polyplast type electrodes: Note limited resistance to chemicals, see «General» below.

Ceramic diaphragm (liquid junction): If protein contamination occurs, the electrode should be immersed in a solution of 0.4% HCl + 5 g/l pepsin for several hours. If blackening of the diaphragm is apparent (silver compounds), the electrode should be immersed in a solution of 0.4% HCl + 76 g/l thiourea.

Sleeve diaphragm: To clean the ground sleeve diaphragm, loosen the sleeve and clean the joint with a brush or similar tool.

Regeneration

pH: Immerse electrode for 10 minutes in 0.1 – 1.0 M NaOH, then 10 minutes in 0.1 – 1.0 M HCl. After regeneration, place it in storage solution for 15 minutes.

ORP: clean metal surfaces with a mildly abrasive medium, e.g. toothpaste or very fine scouring powder.

Disposal

 The design of Hamilton electrodes minimizes environmental impact. According to the EU directive 2002/96/EC the Hamilton sensors should be disposed as waste of electrical and electronic equipment, and not in municipal waste, or it can be sent back to Hamilton for disposal.

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General

The life span of electrodes depends on the customer requirements of response time, zero point and slope. Measurement conditions - especially high temperatures and aggressive samples - may shorten life span. In favorable conditions, life expectancy can be 1 – 3 years at ambient temperature, and 1 – 3 months at 90 °C. Some ageing also occurs during storage.

Polyplast type electrodes have limited chemical resistance. Extended contact with concentrated acids, alkaline solutions, and alcohols should be avoided. Ethers, esters, ketones, aromatic and halogenated hydrocarbons attack the shaft material and should not be allowed to contact the electrode.

Hint: High-quality Hamilton cables extend life time and **yield more stable measurement values!**

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