

# Microlab<sup>®</sup> 600

## Disposable Tip Hand Probe



### Save Time and Money with a Dedicated Diluter

Most tasks in the laboratory involve dilutions of one type or another. A common example is mixing one or more samples and reagents and then bringing the solution up to a final volume with a diluent like water, buffer, etc. Using the Microlab 600 Diluter for these tasks can reduce user-to-user variation, improve consistency, minimize wasted buffer, and save on labor costs.

### When to Use the Disposable Tip Hand Probe (DTHP)

For most analytical applications, the standard Concorde CT Hand Probe provides exceptional performance with no interference from carryover. However, some applications require a higher level of assurance, and that is why the new Disposable Tip Hand Probe was developed. For these sensitive applications, the sample is aspirated into a disposable plastic tip which is thrown away between each sample, eliminating any chance for carryover. Applications that benefit from the DTHP include:

- Forensics—for some applications regulatory considerations make disposable tips the preferred option
- Sterile samples—sterile disposable tips can be used to avoid transferring contamination between sample vessels
- DNA amplification—for applications where a single amplified strand of DNA is enough to impact results

### Benefits of the DTHP:

- Save time and minimize wasted buffer
- Eliminate carryover between samples
- One probe uses 50  $\mu$ L or 1 mL ClickSure<sup>®</sup> tips
- Minimize fatigue with light weight probe and spring-loaded tip ejection
- Scroll wheel quickly changes between favorites
- Configured for left or right-hand use

## How Does it Work?


A tip is installed onto the hand probe.

**Step 1:** The trigger is pressed and the sample is aspirated into the tip using the right syringe while diluent is drawn into the left syringe.

**Step 2:** The trigger is pressed again and the sample is dispensed from the tip.

**Step 3:** The tip is ejected.

**Step 4:** Another trigger dispenses the diluent to complete the dilution.



Additional DTTP Wizards may be available for download at [www.hamiltoncompany.com/microlabwizards](http://www.hamiltoncompany.com/microlabwizards).

### Ordering Information

Description	Part Number
Microlab 600 Advanced Diluter with DTTP	ML625-DTTP*
Disposable Tip Hand Probe (DTTP)	63960-02

### ClickSure® Pipette Tips

Description	Bulk	Racked
50 µL ClickSure Tips, 960 tips/pack	235543	235537
1000 µL ClickSure Tips, 960 tips/pack	235545	235539

\*Ships complete with Disposable Tip Hand Probe, universal valves, fill tubing, Cable Management System, power cord, and the choice of two syringes.

Throughout this document, protected product names may be used without being specifically marked as such. Research use only. Not for use in diagnostics procedures. All rights reserved. All other trademarks are the sole property of their respective owners. © 2024 Hamilton Company. All rights reserved. All trademarks are owned and/or registered by Hamilton Company in the U.S. and/or other countries. Lit. No. L10066 Rev. C – 06/2024

**HAMILTON**

Hamilton Americas & Pacific Rim  
Hamilton Company Inc.  
4970 Energy Way  
Reno, Nevada 89502 USA  
Tel: +1-775-858-3000  
Fax: +1-775-856-7259  
sales@hamiltoncompany.com

Hamilton Europe, Asia & Africa  
Hamilton Central Europe S.R.L.  
str. Hamilton no. 2-4  
307210 Giarmata, Romania  
Tel: +40-356-635-055  
contact.hce.ro@hamilton-ce.com

Web: [www.hamiltoncompany.com](http://www.hamiltoncompany.com)  
USA: 800-648-5950  
Europe: +40-356-635-055

To find a representative in your area, please visit: [www.hamiltoncompany.com/contacts](http://www.hamiltoncompany.com/contacts)

