

Bonaduz, September 30st 2021


Statement: Dry storage time for Incyte-P SU and Conducell-P SU

To whom it may concern,

We would like to confirm that Hamilton tested the storage lifetime of the Incyte-P SU and Conducell-P SU pre-installation for 12 months and post installation for additional 24 months as specified in the datasheet. The test were executed as accelerated test according to ASTM F1980. Based on the test results we can confirm that the shelf life can be divided differently but may not exceed a total of 36 month.

Yours sincerely

Hamilton Bonaduz AG



Peter Pianegonda
Quality Manager Process Analytics



Philipp Arquint
Head of Innovation Process Analytics

Referenced document: 11100197501b (2020-09-17)